Vicks SpeedRead™ Digital Thermometer

Model V912USV5
Model V912BBUSV5
Model V912BBUSV6

This thermometer should be used by adults only. This thermometer is for Household Use only. Do not leave the device unattended with children under the age of 12.

This product contains a lithium button / coin cell battery. If a new or used lithium button/coin cell battery is swallowed or enters the body, it can cause severe internal burns and can lead to death in as little as 2 hours. Always completely secure the battery compartment. If the battery compartment does not close securely, stop using the product, remove the batteries, and keep it away from children. If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.

WARNING: Chemical Burn Hazard. Keep batteries away from children.

READ AND SAVE THESE INSTRUCTIONS

Thank you for purchasing the Vicks SpeedRead™ Digital Thermometer with large screen display and Fever InSight™ feature, an accurate instrument for faster and easier temperature taking. Please read these instructions carefully to help ensure accurate temperatures and safe operation. When used properly, this high quality product will give you years of accurate, dependable service.

INDICATIONS FOR USE / INTENDED USE

The Vicks SpeedRead™ Digital Thermometer is a handheld, battery-powered, predictive, digital stick thermometer intended for the intermittent determination of human body temperature orally, rectally, or under the arm, in a home-use environment for people of all ages (infants, children, and adults).

IMPORTANT SAFETY INFORMATION

- High, prolonged fever requires medical attention especially for young children. Please contact your physician.
- DO NOT** allow children to take their temperatures unattended.
- DO NOT** allow children to walk or run during temperature taking.
- This thermometer is intended for household use only. This product is not intended to diagnose any disease, but is a useful screening tool for temperature. Use of this thermometer is not intended as a substitute for consultation with your physician.
- Consult a physician immediately if temperature is elevated and fever is suspected. High, prolonged fever requires medical attention, especially for young children.
- Parents/guardians should call the physician upon noticing any untoward sign(s) or symptom(s). For example, a child who exhibits irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, even in the absence of fever, or who exhibits a low temperature, may still need to receive medical attention. Children who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.
- Keep out of reach of children. **DO NOT** leave children unattended when obtaining temperature measurement.
- NEVER** use the thermometer for purposes other than those for which it has been intended. Please follow the general safety precautions when using on children.
- DO NOT** use the thermometer if there are signs of damage on the probe or on the thermometer itself. If damaged, **DO NOT** attempt to repair the product.
- This thermometer is not intended for pre-term babies or small-for-gestational age newborn babies (birth to 4 weeks of age). Pre-term is defined as babies born before 37 weeks of gestation. Small-for-gestational-age is defined as a newborn baby, born at 37 weeks or later, with a weight below the 10th percentile for newborn babies of the same gestational age.
- This thermometer is not intended to interpret hypothermic temperatures. If the device displays a temperature of 97.5°F (36.4°C) or less with an LED backlight color of green or white, and the individual is exhibiting atypical symptoms or behaviors, contact your physician or health care professional.
- Temperature elevation may signal a serious illness, especially in neonates and infants, or adults who are old, frail, have a weakened immune system. Please seek professional advice immediately when there is a temperature elevation in: - Neonates and infants under 3 months (Consult your physician immediately if the temperature exceeds 99.4°F or 37.4°C). - Individuals having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer chemotherapy, chronic immunosuppressant treatment, splenectomy). - Individuals who are bedridden (e.g. paraplegia, quadriplegia, stroke, recovering from surgery, nursing home patient, chronic illness). - A transplant recipient (e.g., liver, heart, lung, kidney).
- DO NOT** allow children to walk or run while using the thermometer.
- This thermometer contains small parts that can be swallowed or produce a choking hazard to children.
- ALWAYS** keep the thermometer out of reach of children.

IP67: Dust tight. No ingress of dust. Protected against 30 minute submersion at a depth of 1 meter.

If you have any questions or comments about the Vicks SpeedRead™ Digital Thermometer, call our toll-free Consumer Relations Line: 1-800-477-0457 or visit: VicksThermometers.com/contact-us

CHANGING THE DISPLAY FROM FAHRENHEIT (F) OR CELSIUS (C):

With the thermometer in the OFF position, hold down the center button until the display shows Lo. Release the button. To choose desired scale, press the button again to move between the "F" or "C" icon shown below the Lo. Once the "C" or "F" starts to blink, push the button again to turn off the thermometer. You can now turn on the thermometer and the display will be in Fahrenheit (F) or Celsius (C) based on the scale you set in the above steps. Your thermometer will return to Fahrenheit (F) readings when you replace your battery. You can repeat this process to switch to Celsius (C).

TEMPERATURE READINGS

Temperature readings vary from person to person. Although there is no one "normal" temperature reading, a temperature, taken orally, ranging between 96.4°F and 99.4°F is considered normal. A rectal temperature is generally 1°F higher and an underarm (axillary) temperature is generally 1°F lower than a temperature measurement taken orally. The best method to determine your own normal temperature is to use the thermometer when you are feeling well. Record your reading twice a day (early morning and late afternoon). Take the average of the two temperatures. This is considered your normal body temperature. Any variation from it may indicate some sort of illness and you should consult your physician.

BEFORE USE

Your Vicks SpeedRead™ Digital Thermometer is a 3-in-1 thermometer for oral, rectal, or underarm use. Whichever method is used, it is recommended that for 15 minutes prior to use, avoid eating or drinking any liquids, exercising, taking showers or baths, or smoking.

OPERATING YOUR Vicks SpeedRead™ DIGITAL THERMOMETER

1. TURNING THE UNIT ON

- Turn unit on by pressing the On/Off button (Fig. 1).
- A beep will sound and the display will flash Red/Yellow/Green for 1 second.
- Wait for unit to get ready.
 - The display will perform a self-test (Fig. 2).
 - The display will show last measurement for 2 seconds (Fig. 3).
 - After about 5 seconds, the display will show it's ready (Fig. 4).

Fig. 1

Fig. 2

Fig. 3


Fig. 4

2. SELECTING A TEMPERATURE LOCATION

Oral Use

- YOU MUST COVER THE THERMOMETER'S TIP WITH A SINGLE-USE, DISPOSABLE PROBE COVER.** Probe cover not included. Single-use, disposable probe covers for digital stick thermometers are available for purchase at retail chain stores and pharmacies.
- Place thermometer under the tongue. It is important to place the tip of the thermometer well under the tongue to get a good, accurate reading. (Fig. 5) Place sensor tip in areas marked with "✓". Keep mouth closed and sit still to help ensure an accurate measurement.
- When the peak temperature is reached, 5 beep signals will sound and the "F" or "C" sign will stop flashing. The temperature is now ready to read and will not change when the thermometer is removed.
- Dispose of used probe cover (not included) and wash the tip as recommended.

Fig. 5

Under-the-arm (Axillary) Use

This method is used for babies and young children.

- YOU MUST COVER THE THERMOMETER'S TIP WITH A SINGLE-USE, DISPOSABLE PROBE COVER.** Probe cover not included. Single-use, disposable probe covers for digital stick thermometers are available for purchase at retail chain stores and pharmacies.
- Wipe the underarm with a towel.
- Place the probe tip under the arm so the tip is touching the skin and position your arm next to your body (Fig. 6). With a young child, it is sometimes helpful to hug the child to keep their arm next to their body. This ensures that the room air does not affect the reading.
- When the peak temperature is reached, 5 beep signals will sound and the "F" or "C" sign will stop flashing. The temperature is now ready to read and will not change when the thermometer is removed.
 - Note: This method produces a temperature up to 1°F (0.6°C) lower than a temperature taken orally.**
- Dispose of used probe cover (not included) and wash the tip as recommended.


Fig. 6

Rectal Use

Commonly used for babies, young children, or when it is difficult to take an oral or underarm temperature.

- YOU MUST COVER THE THERMOMETER'S TIP WITH A SINGLE-USE, DISPOSABLE PROBE COVER AND LUBRICATE WITH A WATER SOLUBLE JELLY FOR EASIER INSERTION.** DO NOT use petroleum jelly. Read the label of the lubricant to be sure you are using the right kind. Probe cover not included. Single-use, disposable probe covers for digital stick thermometers are available for purchase at retail chain stores and pharmacies.
- Lay the individual on his/her side. If the individual is an infant, the proper position for a baby is to lay on its stomach with legs hanging down, either across your knees or at the edge of a bed or changing table. This positions the infant's rectum for safe and easy insertion of the thermometer.
- Gently insert the tip of the thermometer NO MORE THAN 1/2 inch into the rectum. If you detect resistance, stop. Hold the thermometer in place during measurement.
- When the peak temperature is reached, 5 beep signals will sound and the "F" or "C" sign will stop flashing. The temperature is now ready to read and will not change when the thermometer is removed.
 - Note: This method produces a temperature up to 1°F (0.6°C) higher than a temperature taken orally.**
- Dispose of used probe cover (not included) and wash and disinfect the tip as recommended.

3. FEVER INSIGHT™ FEATURE

Fever InSight™ helps you to better understand the meaning of your child's temperature with the color indicated on the display. The screen displays Green if no fever, Yellow when the temperature indicates fever and a Red alert for high fever. Consult your doctor when assessing the seriousness of fever.

Color	Range (°F)	Range (°C)	Meaning
Green	96.4-99.4	35.8-37.4	Normal
Yellow	99.5-101.3	37.5-38.5	Fever
Red	101.4+	38.6+	High Fever

4. TURNING THE UNIT OFF

Turn the thermometer off by pressing the On/Off button. If you forget to turn it off, the thermometer will automatically turn off after approximately 1 1/2 minutes.

CARE AND CLEANING

- Clean the probe tip by washing with soap and warm water or by disinfecting with rubbing alcohol (70% isopropyl alcohol).
- Wipe dry with soft cloth.
- DO NOT BOIL OR CLEAN IN DISHWASHER.** This will cause the thermometer to no longer function and will void the warranty.
- DO NOT** wash the unit with any thinner or chemical solvent.

CHANGING THE BATTERIES

The Vicks SpeedRead™ Digital Thermometer features a long life battery. When the "V" appears and flashes in the center top of the display, the battery is exhausted and needs replacing.

- Using a coin, unlock the battery cover and carefully remove the battery with a non-metallic instrument.
- Replace the battery with a 3V, CR1632 silver oxide type or equivalent. Insert the new battery with the positive (+) side facing the battery cover.
- Carefully replace the battery cover and, using a coin, tighten the cover.
- Ensure the battery door is completely secure after replacing the battery.

Remove the battery from the device if it will not be used for an extended period of time in order to avoid damage to the thermometer resulting from a leaking battery.

CARE AND USE

Operating or storing outside these ranges and dropping the thermometer may negatively impact the measurement and potentially damage the thermometer:

- Operating environment must be between 50°F -104°F.
- Operating and storage environment must be between 15-95% relative humidity.
- Store between -13°F and 140°F.
- If core body temperature is lower than ambient temperature the reading may be negatively impacted.

To protect the environment, dispose of product and empty batteries at appropriate collection sites according to national or local regulations. Keep out of reach of young children, elderly and pets.

SPECIFICATIONS

Atmospheric pressure: 70kPa-106kPa
 Battery Life: 500 readings after storage for 1 year or using one time every other day for 2 years.
 Temperature Range: 89.6°F-109.4°F. If less than 89.6°F, "Lo" is displayed. If more than 109.4°F, "Hi" is displayed.
 Beeper: Unit will beep for approximately 8 seconds when peak temperature is reached. For temperatures over 101.3°F, the audible Fever Signal will sound 5 times.
 Error Messages: The display will read "Err" when a malfunction occurs and the unit will power off after 1 minute.
 Accuracy in Water-bath: +/-0.2°F between 96.0°F-107.6°F at room temperature of 71°F.
 Determination of accuracy in a waterbath requires accessing a special test mode. There is no difference or correction between the waterbath reading and the predictive reading for this thermometer. Please contact Kaz for instructions related to use in waterbath mode.
 Conforms to ASTM E1112, Kaz USA, Inc., a Helen of Troy Company, 400 Donald Lynch Boulevard, Suite 300, Marlborough, MA 01752.

1 YEAR LIMITED WARRANTY

You should first read all Instructions before attempting to use this product. Please retain the receipt as proof of and date of purchase. The receipt must be presented when making any claim within the relevant warranty period. Any claim under the warranty will not be valid without a proof of purchase. Your device is warranted for one year (1 year) from date of delivery.

A. This 1 year limited warranty applies to repair or replacement of product found to be defective in material or workmanship. This warranty does not apply to damage resulting from commercial, abusive, unreasonable use or supplemental damage. Defects that are the result of normal wear and tear will not be considered manufacturing defects under this warranty. **KAZ IS NOT LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR THIS PRODUCT IS LIMITED IN DURATION TO THE DURATION OF THIS WARRANTY.**

Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you also may have other rights that vary from jurisdiction to jurisdiction. This warranty applies only to the original purchaser of this product from the original date of delivery.

B. At its option, Kaz will repair or replace this product if it is found to be defective in material or workmanship.

C. This warranty does not cover damage resulting from any unauthorized attempts to repair or to specify a use not in accordance with the instruction manual.

Call us toll-free at 1-800-477-0457 or visit: VicksThermometers.com/contact-us. Please be sure to specify a model number.

NOTE: IF YOU EXPERIENCE A PROBLEM, PLEASE CONTACT CONSUMER RELATIONS FIRST OR SEE YOUR WARRANTY. DO NOT RETURN THE PRODUCT TO THE ORIGINAL PLACE OF PURCHASE. DO NOT ATTEMPT TO OPEN THE THERMOMETER HOUSING YOURSELF, DOING SO MAY VOID YOUR WARRANTY AND CAUSE DAMAGE TO THE PRODUCT OR PERSONAL INJURY.

IP67

IP67 - Dust tight. No ingress of dust. Protected against 30 minute submersion at a depth of 1 meter.

Helen of Troy

© 2026 All rights reserved. Distributed by Kaz USA, Inc., a Helen of Troy Company, 400 Donald Lynch Blvd., Suite 300, Marlborough, MA 01752
www.VicksThermometers.com

Certain trademarks used under license from The Procter & Gamble Company or its affiliates.

Fever InSight™ and SpeedRead™ are trademarks of Kaz Europe Sàrl.

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