
- This thermometer is not intended for pre-term babies or small-for-gestational age newborn babies (birth to 4 weeks of age). Pre-term is defined as babies born before 37 weeks of gestation. Small-for-gestational-age is defined as a newborn baby, born at 37 weeks or later, with a weight below the 10th percentile for newborn babies of the same gestational age.
- This thermometer is not intended to interpret hypothermic temperatures. If the device displays a temperature of 97.5°F (36.4°C) or less with an LED backlight color of green or white, and the individual is exhibiting atypical symptoms or behaviors, contact your physician or health care professional.
- Temperature elevation may signal a serious illness, especially in neonates and infants, or adults who are old, frail, have a weakened immune system. Please seek professional advice immediately when there is a temperature elevation in:
 - Neonates and infants under 3 months (Consult your physician immediately if the temperature exceeds 99.4°F or 37.4°C).
 - Individuals having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer chemotherapy, chronic immunosuppressant treatment, splenectomy).
 - Individuals who are bedridden (e.g. paraplegia, quadriplegia, stroke, recovering from surgery, nursing home patient, chronic illness).
 - A transplant recipient (e.g., liver, heart, lung, kidney).
- **DO NOT** allow children to walk or run while using the thermometer.
- This thermometer contains small parts that can be swallowed or produce a choking hazard to children.
- **ALWAYS** keep the thermometer out of reach of children.

IP22: Protected against solid foreign objects of 12.5 mm diameter and larger. Protected against water drops falling vertically when the device is tilted up to 15°.

CHANGING THE DISPLAY FROM FAHRENHEIT (° F) OR CELSIUS (° C):

1. Make sure the thermometer is turned OFF.
2. Press and hold the blue SITE SELECTION button until C/F scale alternates on the display.
3. Release the SITE SELECTION button when the desired scale is shown (C or F).
4. The device will power off after 2 seconds.
5. After being turned ON again, the most recent temperature stored in memory will be displayed in the newly chosen scale.

TEMPERATURE READINGS

Temperature readings vary from person to person. Although there is no one “normal” temperature reading, a temperature, taken orally, ranging between 97°F and 100°F, is considered normal. A rectal temperature is generally 1°F higher and an underarm (axillary) temperature is generally 1°F lower than a temperature measurement taken orally. The best method to determine your own normal temperature is to use the thermometer when you are feeling well. Record your reading twice a day (early morning and late afternoon). Take the average of the two temperatures. This is considered your normal body temperature. Any variation from it may indicate some sort of illness and you should consult your physician.

BEFORE USE: Your Vicks RapidRead™ Digital Thermometer is a 3-in-1 thermometer for oral, rectal, or underarm use. Whichever method is used, it is recommended that for 15 minutes prior to use, avoid eating or drinking any liquids, exercising, taking showers or baths, or smoking.

CARE AND CLEANING

- The thermometer should be cleaned in between uses, as necessary after oral, underarm, or rectal measurements. Wash entire thermometer with soap and warm water until all visible soil is removed and let air dry. Single-use, disposable probe covers (not included) are not cleanable or reusable and should be discarded after being used for a temperature measurement.
- **DO NOT** boil or clean the thermometer in the dishwasher. Ensure that no liquid enters the interior of the thermometer.
- **NEVER** use abrasive cleaning agents, thinners, or benzene for cleaning, and never soak the device in water or other cleaning liquids.

DISINFECTING

- The thermometer should be disinfected in between uses, as necessary after oral, underarm, or rectal measurements. Wet a paper towel with 10% bleach solution, hold the tip of the device with the wet paper towel for 30 seconds, and then wipe towards the tip for approximately 3 minutes and 30 seconds. Rinse with water for 30 seconds and dry.

CARE & USE

Operating or storing outside these ranges may negatively impact the measurement and potentially damage the thermometer:

- Operating environment must be between 59°F -104°F.
- Operating and storage environment must be between 15-95% relative humidity.
- Store between -13°F and 122°F.
- If core body temperature is lower than ambient temperature the reading may be negatively impacted.

To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations. Keep out of reach of young children, elderly, and pets.

CHANGING THE BATTERIES

The Vicks RapidRead™ Digital Thermometer features a long life battery. When the “**LOW**” appears and flashes in the center of the display, the battery is exhausted and needs replacing.

1. Using a coin, unlock the battery cover and carefully remove the battery with a non-metallic instrument.
2. Replace the battery with a 3V, CR2032 lithium type. Insert the new battery with the positive (+) side facing the battery cover.
3. Carefully replace the battery cover and, using a coin, tighten the cover.

SPECIFICATIONS

Battery Life:	More than 400 measurements or approximately 2 years if used every other day.
Temperature Range:	93.2-109.4°F.
Beeper:	Unit will beep when elevated temperature is reached. For temperatures over 99.4°F oral, 99.2°F underarm, or 100.3°F rectal, the audible Fever Signal will sound 5 times.
Error Messages:	Err will be displayed if too much movement or any odd behavior has taken place during a measurement. The display will read Err when the tip of the thermometer is being moved while taking a reading. The display will read Er1 when a malfunction occurs and the unit will power off after 1 minute. The display will read Er2 when the reading is outside the measurement temperature range. The display will read Er3 when the ambient temperature is outside the operating temperature range. Please see care and use for operating and storage environments.
Accuracy in Water-bath:	+/-0.2°F between 93.2-109.4°F at room temperature of 71°F. Determination of accuracy in a waterbath requires accessing a special test mode. There is no difference or correction between the waterbath reading and the predictive reading for this thermometer. Please contact Kaz for instructions related to use in waterbath mode.

Vicks RapidRead™ Digital Thermometer

VDT972USV4

This thermometer should be used by adults only. This thermometer is for Household Use only. **DO NOT** leave the device unattended with children under the age of 12.

⚠ WARNING: Chemical Burn Hazard. Keep batteries away from children.
This product contains a lithium button / coin cell battery. If a new or used lithium button/coin cell battery is swallowed or enters the body, it can cause severe internal burns and can lead to death in as little as 2 hours. Always completely secure the battery compartment. If the battery compartment does not close securely, stop using the product, remove the batteries, and keep it away from children. If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.

READ AND SAVE THESE INSTRUCTIONS

Thank you for purchasing the Vicks RapidRead™ Digital Thermometer with large screen display and Precision InSight™ feature, a professionally accurate instrument for faster and easier temperature taking. Please read these instructions carefully to ensure accurate temperatures and safe operation. When used properly, this high quality product will give you years of accurate, dependable service.

INDICATIONS FOR USE/INTENDED USE

The Vicks VDT972 RapidRead™ Digital Thermometer is a handheld, battery-powered, predictive, digital stick thermometer intended for the intermittent determination of human body temperature orally, rectally, or under the arm, in a home-use environment for people of all ages (infants, children, and adults).

IMPORTANT SAFETY INFORMATION

High, prolonged fever requires medical attention especially for young children. Please contact your physician. Call or go to Emergency Room immediately if for any of the following:

- Fever in infants under 3 months: Any fever is urgent.
- Fever ≥104°F (40°C) at any age, especially if it repeats.
- Fever accompanied by severe symptoms such as trouble breathing, blue or grey lips/skin, unresponsiveness, seizures, stiff neck, severe headache, severe ear pain, unexplained rash, or repeated vomiting/diarrhea.
- **Regardless of temperature, seek urgent medical attention if the infant is hard to wake, feeding poorly, has breathing difficulties, appears unusually lethargic or inconsolable, or shows signs of dehydration.**

- **DO NOT** allow children to take their temperatures unattended.
- **DO NOT** allow children to walk or run during temperature taking.
- This thermometer is intended for household use only. This product is not intended to diagnose any disease, but is a useful screening tool for temperature. Use of this thermometer is not intended as a substitute for consultation with your physician.
- Consult a physician immediately if temperature is elevated and fever is suspected. High, prolonged fever requires medical attention, especially for young children.
- Parents/guardians should call the physician upon noticing any untoward sign(s) or symptom(s). For example, a child who exhibits irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, even in the absence of fever, or who exhibits a low temperature, may still need to receive medical attention. Children who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.
- Keep out of reach of children. **DO NOT** leave children unattended when obtaining temperature measurement.
- **NEVER** use the thermometer for purposes other than those for which it has been intended. Please follow the general safety precautions when using on children.
- **DO NOT** use the thermometer if there are signs of damage on the probe or on the thermometer itself. If damaged, **DO NOT** attempt to repair the product.

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Fig. 7


Fig. 8

OPERATING YOUR VICKS RAPIDREAD™ DIGITAL THERMOMETER

Tip: For fastest measurements, complete steps 1 and 2 before inserting the thermometer.

1. TURNING THE UNIT ON

- Turn unit on by pressing the On/Off button (Fig. 1).
- A beep will sound and the display will flash Red/Yellow/Green for 1 second.
- Wait for unit to get ready
 - The display will perform a self-test (Fig. 2).
 - The display will show last measurement for 2 seconds (Fig. 3).
 - After last measurement is displayed, the display will enter ready mode (Fig. 4).

2. SELECTING A TEMPERATURE MEASUREMENT SITE

- After the last memory is displayed, the thermometer enters the site selection mode.
- The site selection icons will cycle from oral to underarm to rectal. (Fig. 4).
- When the icon displays the desired measurement site, press and release the blue SITE SELECTION button (Fig. 5).
- The chosen site will flash for 3 seconds and then display the ready screen (Fig. 6). You may now take a reading at the site you selected.
- You may choose an alternate site by pressing the site selection button again (Fig. 5) at any time during the 3 second confirmation period.

WARNING: For accurate readings, you must match the site selected on the thermometer with the measurement site on the patient.

3. TAKING A TEMPERATURE MEASUREMENT

Oral Use:

- An optional single-use, disposable probe cover may be used on the thermometer tip. Disposable probe covers are **not included** and are available for purchase at retail stores and pharmacies.
- Place thermometer under the tongue. It is important to place the tip of the thermometer well under the tongue to get a good, accurate reading (Fig. 7). Place sensor tip in areas marked with “✓”. Keep mouth closed and sit still to help ensure an accurate measurement.
- When the measurement is complete, two beeps will sound. You may then remove the thermometer. The thermometer will then beep with the displayed reading.
- If used, dispose of the probe cover and clean the thermometer tip as recommended.

Under-the-arm (Axillary) Use: This method is used for babies and young children.

- An optional single-use, disposable probe cover may be used on the thermometer tip. Disposable probe covers are **not included** and are available for purchase at retail stores and pharmacies.
- Dry the underarm with a towel before measuring.
- Place the probe tip under the arm so the tip is touching the skin and position your arm next to your body (Fig. 8). With a young child, it is sometimes helpful to hug the child to keep their arm next to their body. This ensures that the room air does not affect the reading.
- When the measurement is complete, two beeps will sound. You may then remove the thermometer. The thermometer will then beep with the displayed reading.
NOTE: This method may produces a temperature 1°F lower than a temperature taken orally.
- If used, dispose of the probe cover and clean the thermometer tip as recommended.

Rectal Use: Commonly used for babies, young children, or when it is difficult to take an oral or underarm temperature.

- An optional single-use, disposable probe cover may be placed over the thermometer tip before use. For easier insertion, lubricate the probe cover or thermometer tip with a water-soluble jelly. Do not use petroleum jelly. Read the lubricant label to ensure you are using the appropriate type. Disposable probe covers are **not included** and are available for purchase at retail stores and pharmacies.
- Lay the individual being measured on his/her side. If the individual is an infant, the proper position for a baby is to lay on its stomach with legs hanging down, either across your knees or at the edge of a bed or changing table. This positions the infant's rectum for safe and easy insertion of the thermometer.
- Gently insert the tip of the thermometer NO MORE THAN 1/2 inch into the rectum. If you detect resistance, stop. Hold the thermometer in place during measurement.
- When the measurement is complete, two beeps will sound. You may then remove the thermometer. The thermometer will then beep with the displayed reading.
NOTE: This method produces a temperature approximately 1°F higher than a temperature taken orally.
- If used, dispose of the probe cover and clean the thermometer tip as recommended.

4. SET A PREFERRED TEMPERATURE MEASUREMENT SITE

To save a preferred temperature measurement site, so that the preferred site will be automatically selected when the thermometer is turned ON, follow these steps.

- Make sure the thermometer is turned OFF.
- Press and hold the blue SITE SELECTION button and, without releasing it, then press and hold the gray POWER button. Release both buttons when the site selection icon appears with 3 flashing dots.
- Press the SITE SELECTION button until the desired temperature measurement site is selected.
- Wait 4 seconds for the site indicator dot to stop flashing.
- Press the POWER button to turn the thermometer OFF.
- The preferred temperature measurement site will be automatically selected when the thermometer is turned ON.
NOTE: With a preferred temperature measurement site selected, an alternate temperature measurement site can still be selected during each use by pressing the blue SITE SELECTION button when the site selection icon appears after the thermometer is turned ON.
- To remove the preferred temperature measurement site setting, open the battery door, remove the battery and reinsert it.
IMPORTANT: To ensure accurate readings, you must always take a temperature reading at the site on the individual that corresponds to the site selected on the thermometer (oral, underarm, or rectal).

5. PRECISION INSIGHT™ FEATURE

Precision InSight™ helps you to better understand the meaning of your child's temperature with the color indicated on the display. The screen displays Green if no fever, Yellow if fever, and Red if high fever. With Precision InSight™, the screen color is customized based upon the measurement site as shown on the chart to the right.

6. TURNING THE UNIT OFF

Turn the thermometer off by pressing the On/Off button. If you forget to turn it off, the thermometer will automatically turn off after approximately 30 seconds.

1 YEAR LIMITED WARRANTY

You should first read all Instructions before attempting to use this product. Please retain the receipt as proof of and date of purchase. The receipt must be presented when making any claim within the relevant warranty period. Any claim under the warranty will not be valid without a proof of purchase. Your device is warrantied for one year (1 year) from date of delivery.

A. This 1 year limited warranty applies to repair or replacement of product found to be defective in material or workmanship. This warranty does not apply to damage resulting from commercial, abusive, unreasonable use or supplemental damage. Defects that are the result of normal wear and tear will not be considered manufacturing defects under this warranty. **KAZ IS NOT LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THIS PRODUCT IS LIMITED IN DURATION TO THE DURATION OF THIS WARRANTY.** Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you also may have other rights which vary from jurisdiction to jurisdiction. This warranty applies only to the original purchaser of this product from the original date of delivery.

B. At its option, Kaz will repair or replace this product if it is found to be defective in material or workmanship.

C. This warranty does not cover damage resulting from any unauthorized attempts to repair or from any use not in accordance with the instruction manual.

Call us toll-free at 1-800-477-0457 or visit: VicksThermometers.com/contact-us. Please be sure to specify a model number.

NOTE: IF YOU EXPERIENCE A PROBLEM, PLEASE CONTACT CONSUMER RELATIONS FIRST OR SEE YOUR WARRANTY. DO NOT RETURN THE PRODUCT TO THE ORIGINAL PLACE OF PURCHASE. DO NOT ATTEMPT TO OPEN THE THERMOMETER HOUSING YOURSELF, DOING SO MAY VOID YOUR WARRANTY AND CAUSE DAMAGE TO THE PRODUCT OR PERSONAL INJURY.

© 2026 All rights reserved.
Distributed by Kaz USA, Inc., a Helen of Troy Company,
400 Donald Lynch Blvd., Suite 300, Marlborough, MA 01752 • www.VicksThermometers.com

Certain trademarks used under license from The Procter & Gamble Company or its affiliates.

Precision InSight™ and RapidRead™ are trademarks owned by Kaz Europe Sàrl.

Protected by U.S. patents D535,898

Made in Vietnam.

Consult your doctor when assessing the severity of fever.

Site	Temperature Ranges and Meaning		
Oral	96.4 - 99.4 °F 35.8 - 37.4 °C	99.5 - 101.3 °F 37.5 - 38.5 °C	101.4+ °F 38.6+ °C
Under-arm	96.2 - 99.2 °F 35.7 - 37.3 °C	99.3 - 101.1 °F 37.4 - 38.4 °C	101.2+ °F 38.5+ °C
Rectal	97.3 - 100.3 °F 36.3 - 37.9 °C	100.4 - 102.2 °F 38.0 - 39.0 °C	102.3+ °F 39.1+ °C
	Normal (Green)	Fever (Yellow)	High Fever (Red)

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Distribuido por Kaz USA, Inc., una empresa de Helen of Troy Company,